

SESAME BM08-XAFS/XRF

Latif Khan (SESAME)

Sergio Lordano (Sirius/LNLS)

Bernd Meyer (Sirius/LNLS)

Simulation Setups

Name	Description
Setup 1	Setup with M2 toroidal without KB mirrors.
Setup 2	Setup with M2 plane and with KB mirrors focusing (we simulated this case here, but I don't have the workspace. Perhaps you can fill these slides with those simulation)
Setup 3A	Latif's workspace sent by email with M2 toroidal creating a secondary source for the KB mirrors, but with large working distance and small demagnification factors.
Setup 3B	Same as before, but modifying the KB mirrors position closer to focus, so that the demagnification factor is larger and the beam is focused down to $\sim 20 \times 20 \mu\text{m}^2$. In this case, the working distance will be smaller ~ 15 cm. This would make better use of the KB mirrors by increasing spatial resolution and flux density at the sample (better signal-to-noise ratio) simultaneously. This can be optimized even further, depending on the beam size that you desire. Note that the KB mirrors specification will also be tighter as the beam focus is decreased, meaning the mirrors will be more expensive.

Setup 1: M2 toroidal

***** SUMMARY OF DISTANCES *****

** DISTANCES FOR ALL O.E. [mm] **

OE	TYPE	p[mm]	q[mm]	src-oe	src-screen
1	EMPTY	6895.0000	0.0000	6895.0000	6895.0000
2	EMPTY	583.0000	0.0000	7478.0000	7478.0000
3	EMPTY	3862.0000	0.0000	11340.0000	11340.0000
4	MIRROR	1338.0000	0.0000	12678.0000	12678.0000
5	CRYSTAL	2500.0000	0.0000	15178.0000	15178.0000
6	CRYSTAL	23.0000	0.0000	15201.0000	15201.0000
7	EMPTY	1500.0000	0.0000	16701.0000	16701.0000
8	MIRROR	1642.0000	9157.0000	18343.0000	27500.0000

** FOCUSING ELEMENTS **

OE	SHAPE	p_foc	q_foc	1/M
4	CYLINDER	?	?	?
8	TOROID	?	?	?

Sum of Alphas [deg]: 180.000000
Sum of Alphas Mod 180 [deg]: 0.000000
Sum of Alphas Mod 360 [deg]: 180.000000

Aperturing Yes ▾

Open slit/Solid stop aperture/slit ▾


Aperture shape Rectangular ▾

Slit width/x-axis [mm]



Slit height/z-axis [mm]

Slit center/x-axis [mm]

Slit center/z-axis [mm]



Setup 1: M2 toroidal



Setup 2: M2 plane and KB mirrors

***** SUMMARY OF DISTANCES *****


** DISTANCES FOR ALL O.E. [mm] **

OE	TYPE	p[mm]	q[mm]	src-oe	src-screen
1	EMPTY	6895.0000	0.0000	6895.0000	6895.0000
2	EMPTY	583.0000	0.0000	7478.0000	7478.0000
3	EMPTY	3862.0000	0.0000	11340.0000	11340.0000
4	MIRROR	1338.0000	0.0000	12678.0000	12678.0000
5	CRYSTAL	2500.0000	0.0000	15178.0000	15178.0000
6	CRYSTAL	23.0000	0.0000	15201.0000	15201.0000
7	EMPTY	1500.0000	0.0000	16701.0000	16701.0000
8	MIRROR	1642.0000	0.0000	18343.0000	18343.0000
9	MIRROR	6957.0000	0.0000	25300.0000	25300.0000
10	MIRROR	1200.0000	1000.0000	26500.0000	27500.0000



** FOCUSING ELEMENTS **

OE	SHAPE	p_foc	q_foc	1/M
4	CYLINDER	?	?	?
9	PARABOLA	10000000000.00	2200.00	4545454.55
10	ELLIPSE	26500.00	1000.00	26.50

Sum of Alphas [deg]: 450.000000
Sum of Alphas Mod 180 [deg]: 90.000000
Sum of Alphas Mod 360 [deg]: 90.000000



Setup 2: M2 plane and KB mirrors



Setup 3A: M2 toroidal + KB mirrors


***** SUMMARY OF DISTANCES *****

** DISTANCES FOR ALL O.E. [mm] **



OE	TYPE	p[mm]	q[mm]	src-oe	src-screen
1	EMPTY	6895.0000	0.0000	6895.0000	6895.0000
2	EMPTY	583.0000	0.0000	7478.0000	7478.0000
3	EMPTY	3862.0000	0.0000	11340.0000	11340.0000
4	MIRROR	1338.0000	0.0000	12678.0000	12678.0000
5	CRYSTAL	2500.0000	0.0000	15178.0000	15178.0000
6	CRYSTAL	23.0000	0.0000	15201.0000	15201.0000
7	EMPTY	1500.0000	0.0000	16701.0000	16701.0000
8	MIRROR	1642.0000	9157.0000	18343.0000	27500.0000
9	EMPTY	0.0000	0.0000	27500.0000	27500.0000
10	MIRROR	2000.0000	0.0000	29500.0000	29500.0000
11	MIRROR	1200.0000	1000.0000	30700.0000	31700.0000

** FOCUSING ELEMENTS **

OE	SHAPE	p_foc	q_foc	1/M
4	CYLINDER	?	?	?
8	TOROID	?	?	?
10	ELLIPSE	2000.00	2200.00	0.91
11	ELLIPSE	3200.00	1000.00	3.20



Setup 3A: M2 toroidal + KB mirrors



Setup 3B: M2 toroidal + KB mirrors

(Stronger focusing – smaller working distance)


***** SUMMARY OF DISTANCES *****

** DISTANCES FOR ALL O.E. [mm] **

OE	TYPE	p[mm]	q[mm]	src-oe	src-screen
1	EMPTY	6895.0000	0.0000	6895.0000	6895.0000
2	EMPTY	583.0000	0.0000	7478.0000	7478.0000
3	EMPTY	3862.0000	0.0000	11340.0000	11340.0000
4	MIRROR	1338.0000	0.0000	12678.0000	12678.0000
5	CRYSTAL	2500.0000	0.0000	15178.0000	15178.0000
6	CRYSTAL	23.0000	0.0000	15201.0000	15201.0000
7	EMPTY	1500.0000	0.0000	16701.0000	16701.0000
8	MIRROR	1642.0000	9157.0000	18343.0000	27500.0000
9	EMPTY	0.0000	0.0000	27500.0000	27500.0000
10	MIRROR	3200.0000	0.0000	30700.0000	30700.0000
11	MIRROR	800.0000	200.0000	31500.0000	31700.0000

** FOCUSING ELEMENTS **


OE	SHAPE	p_foc	q_foc	1/M
4	CYLINDER	?	?	?
8	TOROID	?	?	?
10	ELLIPSE	3200.00	1000.00	3.20
11	ELLIPSE	4000.00	200.00	20.00




Setup 3B: M2 toroidal + KB mirrors

(Stronger focusing – smaller working distance)

KB-VFM





KB-HFM



In each mirror's widget, you can look at the beam footprint on the mirror. This help us choose an adequate mirror dimension. Remember that in Shadow Y is the mirror length (horizontal axis above). The vertical axis is the mirror width. Therefore, VFM can be 100 mm long without losing much flux. Longer HFM could result in larger flux. It's a compromise.

Setup 3B: M2 toroidal + KB mirrors



Setups Comparison

	Setup 1	Setup 2	Setup 3A	Setup 3B
Focus Size FWHM (H x V) [μm ²]	465 x 97	35 x 16	127 x 77	22 x 22
Total Flux [ph/s/100mA]	2.2 x 10 ¹⁰	2.5 x 10 ⁰⁹	2.5 x 10 ¹⁰	2.2 x 10 ¹⁰
Avg. Flux Density [ph/s/μm ² /100mA]	4.9 x 10⁵	4.2 x 10⁶	2.5 x 10⁶	4.5 x 10⁷